


# SeQuent<sup>®</sup> Please

Clinically proven Paclitaxel releasing coronary balloon catheter



Vascular Systems



# SeQuent<sup>®</sup> reLEASE


PACLITAXEL - RELEASING CORONARY BALLOON CATHETER

Restenosis prevention continues to be a challenge in interventional cardiology. SeQuent<sup>®</sup> Please is a novel concept of targeted, polymer-free and homogenous drug delivery designed to overcome the limitations of today's treatment options. Clinically proven, it offers completely new procedural solutions combined with a higher safety profile and superior efficacy.



## SeQuent<sup>®</sup> Please

A new era of treating coronary stenosis



## Matrix Coating Technology


Safe and effective local drug delivery



## The original clinically proven Drug Eluting Balloon

Superior to existing treatment techniques

# A new era of treating coronary stenosis



## The proven choice for your solution oriented angioplasty

### Advantages and benefits

- ▶ Maximum lesion accessibility due to a PTCA balloon-like handling
- ▶ Single shot, short-term Paclitaxel delivery for long-term vessel patency
- ▶ Targeted drug delivery into the vascular wall
- ▶ Bio-absorbable, polymer-free drug coating for the prevention of side effects
- ▶ Homogenous and pattern free drug distribution
- ▶ No stent related edge effects
- ▶ No need for stent-in-stent procedures
- ▶ Unlimited treatment options in case of patients' re-intervention
- ▶ Reduced duration of antiplatelet therapy
- ▶ Single treatment or in combination with a bare metal stent (BMS)

### Beyond the standards

Superior clinical patient outcomes in complex lesions

**ISR - In-stent restenosis**

**SVD - Small vessel disease**

Advanced interventional options for the most challenging cases

**CTO - Chronic Total Occlusion**

**Bifurcations**

**Diabetic patients**

**Long lesions**

## Mode of action




### Pure Paclitaxel coating

A short single dose drug application is only possible if Paclitaxel is bioavailable on the balloon surface. Paclitaxel, if applied directly to the balloon surface, needs a matrix, which allows a reliable release and enables an immediate uptake into the vascular wall. If Paclitaxel is applied as a firm compound the required bioavailability is not given.




### Matrix coating with "spacer"

The SeQuent® Please Matrix Coating is a dispersion of Paclitaxel and lopromide. Lopromide acts as a "spacer" and thereby makes the coating porous and Paclitaxel bioavailable. The hydrophilic character of lopromide and the lipophilic properties of Paclitaxel support the release of the drug from the balloon surface and its delivery into the vascular wall.




Combines targeted drug delivery with PTCA balloon like flexibility and handling.



### Balloon surface in contact with the vascular wall

30-second balloon inflation for effective drug delivery into the vascular wall.



### Paclitaxel migrates into the vascular wall

The matrix dissolves. Paclitaxel migrates into the SMC to prevent proliferation without any residues.


# Matrix Coating Technology



## Safe and effective local drug delivery

### Long-term efficacy with short-term release<sup>1</sup>

After a "single shot" application of Paclitaxel there is a sustained antiproliferative action on SMC over 14 days in absence of cytotoxic effects. Following such a "single shot" drug delivery the Paclitaxel concentration reaches bottom levels in vascular cells after 24 hours.



### Effective and homogenous drug delivery into the vascular wall<sup>2 3 4</sup>

The SeQuent<sup>®</sup> Please drug load is 3 µg/mm<sup>2</sup> balloon surface. A short contact time of only 30 seconds proved to be sufficient to inhibit cell proliferation. Approximately 16% of the total Paclitaxel dose is transferred into the vascular wall.




Stent struts of a DES lead to an inhomogenous patterned drug distribution. About 85 % of the vascular wall is not covered by the struts resulting in low drug tissue level.<sup>2</sup>




Homogenous drug distribution with SeQuent<sup>®</sup> Please technology.<sup>4</sup>

### Reduction of neointimal proliferation<sup>4 5</sup>


The Matrix Coating Technology clearly demonstrated in pre-clinical studies its benefit in suppressing smooth muscle cell (SMC) proliferation. BMS implantation with SeQuent<sup>®</sup> Please leads to a rapid stent endothelialization.<sup>4</sup> The SeQuent<sup>®</sup> Please concept reduces neointimal proliferation up to 4 weeks after intervention (in the porcine model).<sup>5</sup>




Endothelialization at day 5 in the porcine model



Bare Metal Stent (BMS)



SeQuent<sup>®</sup> Please Technology



Sirolimus Eluting Stent

## Evaluation of the Matrix Coating Technology


First in man: ISR I/II study <sup>6 7</sup>

Treatment of coronary in-stent restenosis with a Paclitaxel coated balloon catheter - a prospective, randomized trial

In-stent late loss and restenosis at 6 months and MACE at 2-years:

|                           | Uncoated balloon | Paclitaxel-eluting balloon | p-value |
|---------------------------|------------------|----------------------------|---------|
| In-stent late loss [mm]   | 0.81 ± 0.79      | 0.14 ± 0.46                | 0.001   |
| In-segment late loss [mm] | 0.80 ± 0.79      | 0.11 ± 0.48                | 0.001   |
| In-stent restenosis       | 49 %             | 6 %                        | 0.001   |
| In-segment restenosis     | 51 %             | 6 %                        | 0.001   |
| MACE                      | 46%              | 11%                        | 0.001   |

Major adverse events at 2 years:



### Conclusion:

- ▶ Proven clinical efficacy up to 24 months, safe drug delivery
- ▶ No coating-related adverse events, no late thrombosis
- ▶ Clopidogrel post PCI follow-up treatment for only one month
- ▶ Inhibition of restenosis by drug-coated balloons does not require stent implantation and sustained drug release at the site of injury

- 1 **Axel DI et al.**  
Paclitaxel Inhibits Arterial Smooth Muscle Cell Proliferation and Migration In Vitro and In Vivo Using Local Drug Delivery  
Circulation 1997;96:636-645
- 2 **Hwang CW et al.**  
Physiological Transport Forces Govern Drug Distribution for Stent-Based Delivery,  
Circulation 2001;104:600-605
- 3 **Scheller B, Speck U et al.**  
Paclitaxel Balloon Coating, a Novel Method for Prevention and Therapy of Restenosis  
Circulation 2004;110:810-814
- 4 **Scheller B, Speck U, Böhm M**  
Prevention of restenosis: is angioplasty the answer?  
Heart 2007;93:539-541
- 5 **Speck U, Scheller B et al.**  
Neointima Inhibition: Comparison of Effectiveness of Non-Stent Based Local Drug Delivery and a Drug eluting Stent in Porcine Coronary Arteries  
Radiology 2006;240:411-418
- 6 **Scheller B et al.**  
Treatment of Coronary In-Stent Restenosis with a Paclitaxel-Coated Balloon Catheter  
New England Journal of Medicine 2006;355:2113-24
- 7 **Scheller B, Hehrlein C, Bocksch W, Rutsch W, Haghi D, Dietz U, Böhm M, Speck U**  
Two year follow-up after treatment of coronary in-stent restenosis with a paclitaxel-coated balloon catheter  
Clin Res Cardiol 97:773-781 (2008)
- 8 **Drug-Eluting Balloons May Be Alternative to Drug-Eluting Stents**  
FDAnews Device Daily Bulletin Oct. 26, 2007;Vol. 4;No. 211
- 9 **Unverdorben M et al.**  
Paclitaxel-Coated Balloon Catheter versus Paclitaxel-Coated Stent for the Treatment of Coronary In-stent Restenosis  
Circulation 2009;119:2963-2964

PACLITAXEL - RELEASING CORONARY BALLOON CATHETER



| Balloon diameter | Balloon length | Order number |
|------------------|----------------|--------------|
| 2.5 mm           | 10 mm          | 5022200      |
| 3.0 mm           | 10 mm          | 5022201      |
| 3.5 mm           | 10 mm          | 5022202      |
| 4.0 mm           | 10 mm          | 5022204      |
| 2.0 mm           | 15 mm          | 5022205      |
| 2.5 mm           | 15 mm          | 5022206      |
| 2.75 mm          | 15 mm          | 5022207      |
| 3.0 mm           | 15 mm          | 5022208      |
| 3.5 mm           | 15 mm          | 5022210      |
| 4.0 mm           | 15 mm          | 5022211      |
| 2.0 mm           | 17 mm          | 5022212      |
| 2.5 mm           | 17 mm          | 5022213      |
| 2.75 mm          | 17 mm          | 5022214      |
| 3.0 mm           | 17 mm          | 5022215      |
| 3.5 mm           | 17 mm          | 5022218      |
| 4.0 mm           | 17 mm          | 5022219      |
| 2.0 mm           | 20 mm          | 5022220      |
| 2.5 mm           | 20 mm          | 5022221      |
| 2.75 mm          | 20 mm          | 5022222      |
| 3.0 mm           | 20 mm          | 5022223      |
| 3.5 mm           | 20 mm          | 5022225      |
| 4.0 mm           | 20 mm          | 5022227      |
| 2.5 mm           | 26 mm          | 5022230      |
| 2.75 mm          | 26 mm          | 5022231      |
| 3.0 mm           | 26 mm          | 5022232      |
| 3.5 mm           | 26 mm          | 5022234      |
| 2.5 mm           | 30 mm          | 5022240      |
| 3.0 mm           | 30 mm          | 5022242      |
| 3.5 mm           | 30 mm          | 5022243      |

Order information



# The original clinically proven DEB



## Superior to existing treatment techniques

B. Braun clinical study program: PEPCAD <sup>8</sup> <sup>9</sup>

### PEPCAD I

Treatment of small vessel disease (SVD) in 120 patients. 6-month angiographic data compared to literature data (\*Stone G. Jama 2005; 294:1215-23).

|                           | SeQuent® Please | *PES<br>(Paclitaxel eluting stent) |
|---------------------------|-----------------|------------------------------------|
| In-segment late loss [mm] | 0.16 ± 0.38     | 0.49 ± 0.61                        |
| In-segment restenosis     | 5.5 %           | 31.2 %                             |
| MACE (after 6 months)     | 6.1 %           | 18.9 %                             |
| MACE (after 12 months)    | 6.1 %           |                                    |

#### Conclusion:

These results compare very favourably with previously published results on drug eluting stents (DES) in the treatment of small vessel disease (SVD). SeQuent® Please showed a very low binary restenosis rate (5.5 % vs. 31.2 %) and MACE rate (6.1 % vs. 18.9 %) as compared to Paclitaxel eluting stent data.

### PEPCAD II

131 patients with in-stent restenosis (ISR) were randomized to either the SeQuent® Please or a Paclitaxel eluting stent (PES) treatment group. 6-month angiographic data and the 12-month MACE rate are shown in the following table:

|                           | SeQuent® Please | PES         | p-value |
|---------------------------|-----------------|-------------|---------|
| In-segment late loss [mm] | 0.18 ± 0.41     | 0.39 ± 0.63 | 0.03    |
| In-segment restenosis     | 6.7 %           | 20.4 %      | 0.03    |
| MACE (after 6 months)     | 4.7 %           | 18.3 %      | <0.02   |
| MACE (after 12 months)    | 7.7 %           | 19.0 %      | 0.04    |

#### Conclusion:

Compared to PES, SeQuent® Please was safe and showed clinically superior outcomes after 6 and 12 months.

B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany  
Phone +49 30 689897-0 | Fax +49 30 689897-30 | [www.bbraun.com](http://www.bbraun.com)

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 074 61 95-0 | Fax +49 074 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

 0123

All rights reserved. Technical alterations are possible. This leaflet may be used in no way whatsoever other than for the offering, buying and selling of our products.

No part may be copied or reproduced in any form. In the case of misuse we retain the right to recall our catalogues and pricelists and to take legal action.