

LOW PROFILE DEVICES

VascuFlex[®] 5F

PERIPHERAL SELF-EXPANDING NITINOL STENT

SeQuent[®] Please OTW 018

PERIPHERAL DRUG-COATED BALLOON CATHETER

VascuFlex® 5F


PERIPHERAL SELF-EXPANDING NITINOL STENT

Low Profile Nitinol Stent


- 5F sheath compatibility
- 0.018" guidewire compatibility
- Low chronic outward force (COF)
- High flexibility for optimal vessel apposition

Full range of sizes

- 5 - 10 mm stent diameters
- 20 - 200 mm stent lengths
- 80 cm & 120 shaft lengths
- 180 cm shaft length for transradial/transbrachial approach



Low chronic outward force (COF) for less trauma to the vessel wall



REDUCE THE
MINIMIZE C

LOW PROFILE
0.018" GUIDEW

SeQuent[®] Please OTW

PERIPHERAL DRUG COATED BALLOON CATHETER

LOW PROFILE AND SIMPLE TO USE WITH COMPLICATIONS



Low Profile Drug Coated Balloon

- One system to cover BTK & ATK lesions
- 4F/5F sheath compatibility
- 0.018" guidewire compatibility
- 2.0 – 6.0 mm balloon diameter
- 40 – 120 mm balloon length
- 75 cm & 130 cm shaft length


Resistant DCB coating technology

- Homogeneous coating even in the balloon folds
- Paclitaxel with a drug dose of 3 µg/mm²
- Resveratrol with positive characteristics
- Naturally occurring, plant based substance
- Anti-inflammatory, antioxidant & anti-hypertensive^{2,3}
- Facilitates drug transfer into vessel wall


Proven with strong clinical data over 24 months

- 153 patients (78 DCB vs. 75 POBA)
- Challenging real world lesions (23.5 % TASC C/D lesions)
- Mean lesion length 13.2 ± 10.4 cm
- Predilatation in only 55.6 % of cases


6 month late lumen loss (LLL)¹



24 month TLR rate



24 month mean walking distance¹



LOW PROFILE DEVICES FOR 4F/5F SHEATH

¹ Albrecht T et al. Cardiovasc Intervent Radiol (2018). doi:10.1007/s00270-018-1940-1

² Ciccone G. et al. in Curr Med Chem. 2013 ;20(10):1323-31.

³ Magyar K. et al. Clin Hemorheol Microcirc. 2012; 50(3):179-87.

VascuFlex® 5 F

		5 mm	6 mm	7 mm	8 mm	9 mm	10 mm
VascuFlex® 5 F 80 cm Shaft Length	20 mm	5505060	5505070	5505080	5505090	5505100	5505110
	40 mm	5505061	5505071	5505081	5505091	5505101	5505111
	60 mm	5505062	5505072	5505082	5505092	5505102	5505112
	80 mm	5505063	5505073	5505083	5505093	5505103	5505113
	100 mm	5505064	5505074	5505084	5505094	5505104	5505114
	120 mm	5505065*	5505075*	5505085*	5505095*	5505105*	5505115*
	150 mm	5505066*	5505076*	5505086*	5505096*		
	180 mm	5505067*	5505077*	5505087*			
	200 mm	5505068*	5505078*	5505088*			

VascuFlex® 5 F 120 cm Shaft Length	20 mm	5505000	5505010	5505020			
	40 mm	5505001	5505011	5505021	5505031	5505041	5505051
	60 mm	5505002	5505012	5505022	5505032	5505042	5505052
	80 mm	5505003	5505013	5505023	5505033	5505043	5505053
	100 mm	5505004	5505014	5505024	5505034	5505044	5505054

VascuFlex® 5F TR/TB 180 cm Shaft Length	40 mm		5505121	5505124	5505127		
	60 mm		5505122	5505125	5505128		
	80 mm		5505123	5505126	5505129		

*available on request

SeQuent® Please OTW 018

Shaft length 75 cm	Shaft length 130 cm	Balloon diameter (mm)	Balloon length (mm)	Sheath compatibility (F)	Nominal pressure (atm)	Rated burst pressure (atm)
	18220040	2.0 mm	40 mm	4 F	6 atm	16 atm
	18220080	2.0 mm	80 mm	4 F	6 atm	15 atm
	18220120	2.0 mm	120 mm	4 F	6 atm	15 atm
	18225040	2.5 mm	40 mm	4 F	6 atm	16 atm
	18225080	2.5 mm	80 mm	4 F	6 atm	15 atm
	18225120	2.5 mm	120 mm	4 F	6 atm	15 atm
	18230040	3.0 mm	40 mm	4 F	6 atm	16 atm
	18230080	3.0 mm	80 mm	4 F	6 atm	15 atm
	18230120	3.0 mm	120 mm	4 F	6 atm	15 atm
18135040	18235040	3.5 mm	40 mm	4 F	6 atm	16 atm
18135080	18235080	3.5 mm	80 mm	4 F	6 atm	15 atm
18135120	18235120	3.5 mm	120 mm	4 F	6 atm	15 atm
18140040	18240040	4.0 mm	40 mm	4 F	6 atm	16 atm
18140060	18240060	4.0 mm	60 mm	4 F	6 atm	16 atm
18140080	18240080	4.0 mm	80 mm	4 F	6 atm	15 atm
18140120	18240120	4.0 mm	120 mm	4 F	6 atm	15 atm
18150040	18250040	5.0 mm	40 mm	4 F	6 atm	14 atm
18150060	18250060	5.0 mm	60 mm	4 F	6 atm	14 atm
18150080	18250080	5.0 mm	80 mm	4 F	6 atm	13 atm
18150120	18250120	5.0 mm	120 mm	4 F	6 atm	13 atm
18160040	18260040	6.0 mm	40 mm	5 F	6 atm	14 atm
18160060	18260060	6.0 mm	60 mm	5 F	6 atm	14 atm
18160080	18260080	6.0 mm	80 mm	5 F	6 atm	12 atm
18160120	18260120	6.0 mm	120 mm	5 F	6 atm	12 atm

Distributor

B. Braun Melsungen AG | Sieversufer 8 | 12359 Berlin | Germany
Phone +49 30 568207-300 | Fax +49 30 568207-130 | www.bbraun.com

Manufacturer acc. to MDD 93/42/EEC of SeQuent® Please OTW and VascuFlex® Multi-LOC are B. Braun Melsungen AG, Carl-Braun-Str. 1, 34212 Melsungen, Germany. The product trademarks 'SeQuent®' and 'VascuFlex®' are registered trademarks of B. Braun Melsungen AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.