



AESCULAP® Power Systems

Catalog 2024

ELAN 4 and Acculan 4

Content

ELAN 4

- › ELAN 4 electro
- › ELAN 4 Accessories

Acculan 4

- › Acculan 4
- › Acculan 4 Small Drill
- › Acculan 4 Electro
- › Acculan 4 Accessories

Processing

- › Sterilit® Power Systems
- › Processing

Others

- › Tool Attachments Overview
 - › Signs & Symbols
 - › Brochures
 - › Index
-

All measurements and technical data listed in this catalogue are approximate as they are subject to normal manufacturing tolerances and to technical modifications.

AESCLAP® ELAN 4

Precision. Performance. Pleasure.

With its technical concept ELAN 4 addresses important requirements of the surgeon and all parties involved in handling power systems in the hospital.

- Convenient and reliable grip
- Powerful, precise and reliable work
- Easy to set up and simple to use



AESCULAP® ELAN 4

ELAN 4 electro

Electric Highspeed System



ELAN 4 electro

Basic components

GA800

ELAN 4 electro control unit

with irrigation pump; touch screen;
pump flow 0–65 ml/min

100–240 V AC; 50/60 Hz;
380 x 330 x 200 mm; 9,500 g

consisting parts:
GD412804 bottle holder
power cord not included



GA808

ELAN 4 electro foot control

protection type IPX8; anaesthetic
proof class AP

5.0 m;
245 x 144 x 76 mm; 2,240 g

consisting parts:
GA808207 rail



GA810

ELAN 4 electro wireless foot control

protection type IPX8; anaesthetic
proof class AP

271 x 266 x 85 mm; 3,800 g

Requires control unit with software
version V3.00.

consisting parts:
GA810202 rail
3 batteries Type D/LR20/Mono



GA806

GB073R

ELAN 4 electro motor cable

for use with ELAN 4 electro applied
parts

4.0 m; 270 g

consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



AESCULAP® ELAN 4

ELAN 4 electro

Basic components, perforator driver & lowspeed motor

GA804

GB073R

ELAN 4 electro motor cable with hand switch
for use with ELAN 4 electro applied parts
4.0 m; 374 g
consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



GA805

GB073R

ELAN 4 electro motor cable with hand control
for use with ELAN 4 electro applied parts
4.0 m; 369 g
consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



GA822

GB076R

ELAN 4 electro perforator driver
for use with cranial perforators with Hudson shank
180 W; 0 – 1,200 rpm; 250 Ncm; 166 x 35 mm; 410 g




GA824

GB077R

ELAN 4 electro lowspeed motor with intra coupling
180 W; 0 – 40,000* rpm; 10 Ncm; 123 x 23 mm; 150 g

* Requires control unit with software version V3.00. With lower versions the motor runs max. 20,000 rpm.



ELAN 4 electro

Craniotome & attachments

GA849

GB085R


ELAN 4 electro craniotome and multifunction handpiece
2-ring coding; for use with dura guards and holding sleeves
GB941R – GB947R; spray nozzle GA259SU (not included)
140 W; 0 – 80,000 rpm; 2.2 Ncm; 156 x 17 mm; 121 g



GB941R

GB719R

ELAN 4 dura guard pediatric
fixed; 1 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB942R

GB719R

ELAN 4 dura guard standard
fixed; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB943R

GB719R

ELAN 4 dura guard long
fixed; 3 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB945R

GB719R

ELAN 4 holding sleeve
for craniotome and multifunction handpieces GA749, GA849



GB947R





GB719R

ELAN 4 dura guard standard
turnable; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849








AESCULAP® ELAN 4

ELAN 4 electro 1-ring highspeed handpieces

<div>GA861</div> <div>GB083R</div> <div>ELAN 4 electro 1-ring handpiece L4</div> <div>1-ring coding; spray nozzle GA261SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 152 x 16 mm; 80 g</div>	 <div>~ 4 cm</div>
<div>GA862</div> <div>GB083R</div> <div>ELAN 4 electro 1-ring handpiece L7</div> <div>1-ring coding; spray nozzle GA262SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 177 x 16 mm; 83 g</div>	 <div>~ 7 cm</div>
<div>GA863</div> <div>GB083R</div> <div>ELAN 4 electro 1-ring handpiece L10</div> <div>1-ring coding; spray nozzle GA263SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 206 x 16 mm; 87 g</div>	 <div>~ 10 cm</div>
<div>GA864</div> <div>GB083R</div> <div>ELAN 4 electro 1-ring handpiece L13</div> <div>1-ring coding; spray nozzle GA264SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 237 x 16 mm; 91 g</div>	 <div>~ 13 cm</div>

ELAN 4 electro 2-ring highspeed handpieces

<div>GA865</div> <div>GB083R</div> <div>ELAN 4 electro 2-ring handpiece L7</div> <div>2-ring coding; spray nozzle GA265SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 176 x 16 mm; 89 g</div>	 <div>~ 7 cm</div>
<div>GA866</div> <div>GB083R</div> <div>ELAN 4 electro 2-ring handpiece L10</div> <div>2-ring coding; spray nozzle GA266SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 206 x 16 mm; 95 g</div>	 <div>~ 10 cm</div>
<div>GA867</div> <div>GB083R</div> <div>ELAN 4 electro 2-ring handpiece L13</div> <div>2-ring coding; spray nozzle GA267SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 236 x 16 mm; 100 g</div>	 <div>~ 13 cm</div>
<div>GA868</div> <div>GB083R</div> <div>ELAN 4 electro 2-ring handpiece L17</div> <div>2-ring coding; spray nozzle GA268SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 276 x 16 mm; 108 g</div>	 <div>~ 17 cm</div>
<div>GA869</div> <div>GB083R</div> <div>ELAN 4 electro 2-ring handpiece L22</div> <div>2-ring coding; spray nozzle GA269SU (not included) 140 W; 0-80,000 rpm; 2.2 Ncm; 326 x 16 mm; 118 g</div>	 <div>~ 22 cm</div>

AESCULAP® ELAN 4

ELAN 4 electro
MIS handpiece & handpiece shafts



GA860 **GB084R**

ELAN 4 electro MIS handpiece
for use with MIS handpiece shafts
GB920R - GB927R
140 W; 0-80,000 rpm;
2.2 Ncm; 118 x 17 mm; 100 g
consisting parts:
GB600870 Sterilit® spray adapter for
ELAN 4 MIS handpiece shafts



GB920R **GB723R**

**ELAN 4 MIS handpiece shaft
L10 straight**
spray nozzle GB796SU (not included)



GB921R **GB723R**

**ELAN 4 MIS handpiece shaft
L10 curved**
spray nozzle GB796SU (not included)



GB922R **GB723R**

**ELAN 4 MIS handpiece shaft
L10 heavy curved**
spray nozzle GB796SU (not included)



GB925R **GB723R**

**ELAN 4 MIS handpiece shaft
L13 straight**
spray nozzle GB797SU (not included)



GB926R **GB723R**

**ELAN 4 MIS handpiece shaft
L13 curved**
spray nozzle GB797SU (not included)



GB927R **GB723R**

**ELAN 4 MIS handpiece shaft
L13 heavy curved**
spray nozzle GB797SU (not included)



AESCULAP® ELAN 4


ELAN 4 electro Saws



GA836

GB086R

ELAN 4 electro micro sagittal saw
spray nozzle GA258SU (not included)
180 W; 0–20,000 rpm; 161 x 23 mm;
∠ 7°; 190 g



GA837

GB087R

ELAN 4 electro micro reciprocating saw
spray nozzle GA258SU (not included)
180 W; 0–20,000 rpm; 184 x 23 mm;
∅ 1.3 mm; 198 g



GA831

GB088R

ELAN 4 electro sagittal saw
spray nozzle GA258SU (not included)
180 W; 0–20,000 rpm; 129 x 23 mm;
∠ 4.5°; 220 g



GA832

GB087R

ELAN 4 electro reciprocating saw
spray nozzle GA258SU (not included)
180 W; 0–20,000 rpm; 194 x 23 mm;
∅ 2.5 mm; 238 g



GA833

GB088R

ELAN 4 electro transversal saw
spray nozzle GA258SU (not included)
180 W; ∠ 10°; 0–20,000 rpm; 10 Ncm;
130 x 23 mm; ∠ 10°; 185 g
to be used with saw blade shafts
GP580R, GP582R, GP586R and GP590R



AESCULAP® ELAN 4

ELAN 4 electro
Drill



GA844

GB482R

ELAN 4 electro drill

with fixed motor cable, oscillating and tapping mode

200 W; 3.2 mm;
122 x 130 x 27 mm; 800 g

consisting parts:

GA344244 K-wire protection guide

TA011944 cleaning brush

GA344211 rinsing adapter

GB600880 Sterilit® spray adapter

Requires control unit with software version V3.00.



GB870R

GB716R

Intramedullary reamer attachment

with Jacobs chuck 0.5 - 7.4 mm,
not for use with acetabulum reamers

0 - 320 rpm; 6.5 Nm; 3.2 mm;
105 x 29 mm; 230 g

consisting parts:

key GA031R



GB871R

GB716R

Intramedullary reamer attachment

with large AO chuck,
not for use with acetabulum reamers

0 - 320 rpm; 6.5 Nm; 3.2 mm;
88 x 29 mm; 170 g



GB872R

GB716R

Intramedullary reamer attachment

with Hudson/Zimmer chuck,
not for use with acetabulum reamers

0 - 320 rpm; 6.5 Nm; 3.2 mm;
86 x 25 mm; 150 g



AESCULAP® ELAN 4

ELAN 4 electro Attachments

GB881R GB716R

Drill attachment
with Jacobs chuck 0.5 – 7.4 mm
0 – 1,250 rpm; 2 Nm; 3.2 mm;
89 x 29 mm; 163 g
consisting parts:
key GA031R



GB882R GB716R

Drill attachment
with Jacobs chuck 0.5 – 4 mm
0 – 1,250 rpm; 2 Nm; 3.2 mm;
76 x 22 mm; 100 g
consisting parts:
key GA062R



GB884R GB716R

Drill attachment
with small AO chuck
0 – 1,250 rpm; 2 Nm; 2.8 mm;
60 x 22 mm; 66 g



GB886R GB716R

Drill attachment
with AESCULAP® hexagonal chuck
0 – 1,250 rpm; 2 Nm; 3.2 mm;
57 x 24 mm; 67 g



GB887R GB716R

Drill attachment
with Dental shank
0 – 1,250 rpm; 2 Nm; 1.8 mm;
56 x 22 mm; 50 g



GB878R GB716R

Drill attachment
for Synthes radiolucent drive 511.300
0 – 1,250 rpm; 2 Nm;
111 x 27 mm; 181 g



GB894R GB716R

K-wire attachment
for K-wires from 0.6 – 3.2 mm
0 – 1,250 rpm; 3.2 mm;
92 x 27 x 124 mm; 179 g



GB891R GB716R

Sagittal saw attachment
for "L-sagittal" sawblades
0 – 15,750 rpm; 74 x 33.5 mm;
4.5°; 179 g



GB892R GB716R

Reciprocating saw attachment
for reciprocating sawblades
0 – 15,750 rpm; 135 x 24.5 mm;
2.5 mm; 166 g



AESCULAP® ELAN 4



ELAN 4 electro Attachments

GB896R

GB716R

Screw attachment

with small AO chuck

0 - 165 rpm; 12 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  235 g





GB897R

GB716R

Screw attachment






with AESCULAP® hexagonal chuck







0 - 165 rpm; 12 Nm;  1.7 mm;
92 x 28 x 129.5 mm;  231 g



AESCULAP® ELAN 4

ELAN 4 accessories

GB861P GB727R	
GB863P GB727R	
GB862P GB727R	
GB866R	
TE780	
Power cord Europe, 1.5 m	
TE730	
Power cord Europe, 5.0 m	
TE734	
Power cord UK, 5.0 m	
TE735	
Power cord USA, Canada, Japan, 3.5 m	

GA395SU sterile  	
ELAN 4 SPRAY NOZZELS sterile 	
GA259SU for craniotome attachments	
GA261SU for 1-ring handpiece L4	
GA262SU for 1-ring handpiece L7	
GA263SU for 1-ring handpiece L10	
GA264SU for 1-ring handpiece L13	
GA265SU for 2-ring handpiece L7	
GA266SU for 2-ring handpiece L10	
GA267SU for 2-ring handpiece L13	
GA268SU for 2-ring handpiece L17	
GA269SU for 2-ring handpiece L22	
GB796SU for MIS handpiece shafts L10	
GB797SU for MIS handpiece shafts L13	
GA258SU for saws	

AESCULAP® Acculan 4

Look. Feel. Experience.

Acculan 4 the battery-operated power system for orthopaedics, traumatology, cardiothoracic surgery and dermatology

- Powerful
- Longevity – Exclusion of battery and electronic control from sterilization process
- Easy handling – Titanium housing allows mechanical alkaline reprocessing



AESCULAP® Acculan 4

Drill & attachments

GA330 GB495R

Acculan 4 drill and reamer
250 W; 4.0 mm;
207 x 121 x 55 mm; 1,250 g
consisting parts:
GA678 sterile funnel
GA672207 K-wire protection sleeve
TA011944 cleaning brush
GA672211 rinsing adapter
GB600840 Sterilit® spray adapter
Acculan

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GB620R GB496R

Acculan 4 drill attachment
with Jacobs chuck 0.5 – 7.4 mm
0 – 1,000 rpm; 5.0 Nm; 4.0 mm;
133 x 35 mm; 397 g
consisting parts:
key GA031R



GB621R GB496R

Acculan 4 drill attachment
with Jacobs chuck 0.5 – 7.4 mm
0 – 1,000 rpm; 5.0 Nm; 4.0 mm;
118 x 35 mm; 300 g
consisting parts:
key GA031R



GB623R GB496R

Acculan 4 drill attachment
with small AO chuck
0 – 1,250 rpm; 4.0 Nm; 2.8 mm;
90 x 35 mm; 210 g



Attachments

GB639R GB496R

Acculan 4 drill attachment
with Jacobs chuck 0.5 – 4.0 mm
0 – 1,250 rpm; 4.0 Nm; 4.0 mm;
105 x 35 mm; 243 g
consisting parts:
key GA062R



GB645R GB496R

Acculan drill attachment
for Synthes radiolucent drive 511.300
0 – 900 rpm; 110 x 35 mm; 225 g



GB663R GB496R

Acculan 4 drill attachment
with trinkle chuck
0 – 1,000 rpm; 5.0 Nm; 4.0 mm;
93 x 35 mm; 217 g



GB664R GB496R

Acculan 4 drill attachment
with AESCULAP® hexagonal chuck
0 – 1,000 rpm; 5.0 Nm; 4.0 mm;
86 x 35 mm; 225 g



GB665R GB496R

Acculan 4 drill attachment
with Hudson / Zimmer chuck
0 – 1,000 rpm; 5.0 Nm; 4.0 mm;
100 x 35 mm; 230 g



AESCULAP® Acculan 4

Attachments

GB657R

GB496R

Acculan 4 intramedullar attachment

with large AO chuck

0–370 rpm; 12.0 Nm; 4.0 mm;
108 x 35 mm; 300 g



GB667R


GB496R

Acculan 4 reamer attachment

with Jacobs chuck 0.5–7.4 mm chuck

0–250 rpm; 19.0 Nm; 4.0 mm;
133 x 35 mm; 280 g

consisting parts:
key GA031R



GB668R

GB496R

Acculan 4 reamer attachment

with large AO chuck

0–250 rpm; 19.0 Nm; 4.0 mm;
108 x 35 mm; 300 g



GB669R

GB496R

Acculan 4 reamer attachment

with Hudson/Zimmer chuck

0–250 rpm; 19.0 Nm; 4.0 mm;
118 x 35 mm; 300 g



GB670R

GB496R

Acculan 4 reamer attachment

with Harris chuck

0–250 rpm; 19.0 Nm; 4.0 mm;
106 x 35 mm; 290 g



GB660R

GB496R

Acculan sagittal saw attachment

for Acculan rapid action saw blades
with a maximum useable length of
50 mm

0–17,000 rpm; 4°47';
89 x 36 x 44 mm; 270 g



GB641R

GB496R

Acculan K-wire attachment

for K-wires from 0.6–4.0 mm

0–1,250 rpm; 93 x 36 x 137 mm;
285 g



AESCULAP® Acculan 4

Oscillating saws

GA331 GB495R

Acculan 4 oscillating saw
250 W; ↔ 0 – 13,000 rpm; ∠ 4.1°;
211 x 177 x 55 mm; 1,600 g
consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA336 GB498R

Acculan 4 oscillating saw straight
250 W; ↔ 0 – 13,000 rpm; ∠ 4.1°;
290 x 55 x 84 mm; 1,480 g
consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA337 GB498R

Acculan 4 oscillating sternum saw
250 W; ↔ 0 – 13,000 rpm; ∠ 4.1°;
290 x 55 x 84 mm; 1,480 g
consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



Reciprocating saw

GA334 GB495R

Acculan 4 reciprocating saw
250 W; ↔ 0 – 15,000 rpm; ↑ 3.2 mm;
207 x 180 x 55 mm; 1,540 g
consisting parts:
GA678 sterile funnel
GB440215 cleaning plate
GA674890 cleaning adapter

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GB436R GB686R

Protection guide
for use with GA334 and saw blade
GC761R



AESCULAP® Acculan 4

Dermatome

GA340

GB498R

Acculan 4 dermatome

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.2 - 1.2 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; \leftrightarrow 0 - 6,500 rpm; \updownarrow 3.1 mm;
277 x 113 x 83 mm; \blacksquare 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



GA341

GB498R

Acculan 4 dermatome 0.1 mm

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.1 - 1.1 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; \leftrightarrow 0 - 6,500 rpm; \updownarrow 3.1 mm;
277 x 113 x 83 mm; \blacksquare 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



BA720R

GB688R

Skin graft mesher

for use in combination with skin graft
carrier plates BA721, BA722, BA723
(factors 1.5, 3.0, 6.0)

150 x 195 x 105 mm; \blacksquare 4,800 g

consisting parts:
BA725R cutting roller
BA726R ratchet wrench



AESCULAP® Acculan 4

Look. Feel. Experience.

Acculan 4 small drill – the smallest member of the Acculan 4 family can be used for small fragment Osteosynthesis as well as for reamed intramedullary nailing.



AESCULAP® Acculan 4

Small drill
Drill & attachments

GA344 GB482R

Acculan 4 small drill
with oscillating and tapping mode;
not for use with acetabulum reamers
85 W; 3.2 mm;
177 x 124 x 55 mm; 840 g
consisting parts:
GA678 sterile funnel
GA344244 K-wire protection guide
TA011944 cleaning brush
GA344211 rinsing adapter
GB600880 Sterilit® spray adapter

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.



Measurements and weight are including battery and lid.

GB870R GB716R

Intramedullary reamer attachment
with Jacobs chuck 0.5 - 7.4 mm,
not for use with acetabulum reamers
0 - 320 rpm; 5 Nm; 3.2 mm;
105 x 29 mm; 230 g
consisting parts:
key GA031R



GB871R GB716R

Intramedullary reamer attachment
with large AO chuck,
not for use with acetabulum reamers
0 - 320 rpm; 5 Nm; 3.2 mm;
88 x 29 mm; 170 g



GB872R GB716R

Intramedullary reamer attachment
with Hudson/Zimmer chuck,
not for use with acetabulum reamers
0 - 320 rpm; 5 Nm; 3.2 mm;
86 x 25 mm; 150 g



Small drill
Attachments

GB881R GB716R

Drill attachment
with Jacobs chuck 0.5 - 7.4 mm
0 - 1,250 rpm; 1.5 Nm; 3.2 mm;
89 x 29 mm; 163 g
consisting parts:
key GA031R



GB882R GB716R

Drill attachment
with Jacobs chuck 0.5 - 4 mm
0 - 1,250 rpm; 1.5 Nm; 3.2 mm;
76 x 22 mm; 100 g
consisting parts:
key GA062R



GB884R GB716R

Drill attachment
with small AO chuck
0 - 1,250 rpm; 1.5 Nm; 2.8 mm;
60 x 22 mm; 66 g



GB886R GB716R

Drill attachment
with AESCULAP® hexagonal chuck
0 - 1,250 rpm; 1.5 Nm; 3.2 mm;
57 x 24 mm; 67 g



GB887R GB716R

Drill attachment
with Dental shank
0 - 1,250 rpm; 1.5 Nm; 1.8 mm;
56 x 22 mm; 50 g



AESCULAP® Acculan 4

Small drill Attachments

GB878R

GB716R

Drill attachment

for Synthes radiolucent drive 511.300

0 – 1,250 rpm; 1.5 Nm;
111 x 27 mm; 181 g



GB894R

GB716R

K-wire attachment

for K-wires from 0.6 – 3.2 mm

0 – 1,250 rpm; 3.2 mm;
92 x 27 x 124 mm; 179 g



GB891R

GB716R

Sagittal saw attachment

for "L-sagittal" sawblades

0 – 15,750 rpm; 74 x 33.5 mm;
4.5°; 179 g



GB892R

GB716R

Reciprocating saw attachment

for reciprocating saw blades

0 – 15,750 rpm; 135 x 24.5 mm;
166 g; 2.5 mm



GB896R

GB716R

Screw attachment

with small AO chuck

0 – 165 rpm; 10 Nm; 1.7 mm;
95 x 28 x 129.5 mm; 235 g



GB897R

GB716R

Screw attachment

with AESCULAP® hexagonal chuck

0 – 165 rpm; 10 Nm; 1.7 mm;
92 x 28 x 129.5 mm; 231 g




AESCULAP® Acculan 4

Charger & batteries

<div>GA320</div> <div>Acculan 4 charger</div> <div>for Acculan Li-Ion and NiMH batteries GA346, GA676; 4 charging bays; automatic battery diagnosis function 100 - 240 V AC; 50 / 60 Hz; 410 x 120 x 160 mm; 2,300 g power cord not included</div> <div></div>		
<div>GA346</div> <div>Acculan 4 Li-Ion battery short</div> <div>only for GA340, GA341 and GA344; with integrated electronic motor control 9.9 V; 10.89 Wh 120 x 43 x 50 mm; 265 g</div> <div></div>	<div>GA676</div> <div>Acculan NiMH battery</div> <div>only for GA330, GA331, GA334, GA336, GA337; with integrated electronic motor control 9.6 V; 1.95 Ah 140 x 43 x 50 mm; 425 g</div> <div></div>	<div>GA679</div> <div>Acculan battery removal device</div> <div></div>
<div>GA675</div> <div>Acculan lid</div> <div>90 g</div> <div></div>	<div>GA678</div> <div>Acculan sterile funnel</div> <div></div>	<div>GB487R</div> <div>GB488R</div>

AESCULAP® Acculan 4 Electro


Look. Feel. Experience.



Console & batteries

GA877

Acculan Electro console
power cord not included
100 - 240 V AC; 50 / 60 Hz;
255 x 217 x 90 mm; 2,930 g



GA878

GB486R

Acculan Electro connection cable
for use with Acculan 4 handpieces
3.0 m; 410 g



GA875

GB485R

Acculan Electro lid



GA876

Acculan Electro adapter
120 x 49 x 42 mm; 165 g



GA678



GB488R

Acculan sterile funnel



AESCULAP® Acculan 4 Accessories

Adapters, keys & power cord

GB176R GB690R Zimmer adapter tool: Zimmer motor: Hudson/Zimmer only in combination with GB630R, GB656R, GB665R, GB669R or GB872R 3.0 mm	GB414R JG344 AESCULAP® hex adapter tool: AESCULAP® hexagonal motor: triangular 3.5 mm	GB628R JG344 Synthes DHS/DCS adapter tool: Synthes DHS/DCS motor: AO large 4.0 mm
		
GB184R JG344 Harris adapter tool: Harris motor: AO large 3.4 mm	XF457R JG344 DIN adapter tool: DIN motor: AO large 3.7 mm	GA031R JG342 + JG344 Key for large Jacobs chuck suitable for GB620R, GB621R, GB667R, GB870R, GB881R
		
GA062R JG342 Key for small 4 mm Jacobs chuck suitable for GB639R, GB882R		
		
TE780 Power cord Europe, 1.5 m	TE734 Power cord UK, 5.0 m	
TE730 Power cord Europe, 5.0 m	TE735 Power cord USA, Canada, Japan, 3.5 m	

AESCULAP® Processing

Sterilit® power systems



AESCULAP® Processing

Sterilit® power systems

GB600

6

Sterilit® power systems

oil spray; 300 ml



GB600860

Sterilit® power systems spray adapter

for ELAN 4 electro



GB600870

Sterilit® power systems spray adapter

for ELAN 4 MIS handpiece shafts



GB600840

Sterilit® power systems spray adapter

for Acculan 4



GB600880

Sterilit® power systems spray adapter

for Acculan 4 small drill and ELAN 4 electro drill



AESCULAP® Mechanical Processing

- The solution for storage and mechanical reprocessing of our Aesculap® Power Systems.
- For detailed information please refer to the instructions for use.



AESCULAP® Mechanical Processing


ECCOS® sets

ELAN 4

GB261R

ELAN 4 electro ECCOS® set

56 mm
consisting parts:
GB073R
GB076R
GB083R
GB085R
GB719R
GB072R
JF212R



GB262R

ELAN 4 electro ECCOS® set for drill GA844

106 mm
consisting parts:
GB072R
GB073R
GB482R
GB715R
GB716R
2x GB690R
JG342
JG344



GB073R

ECCOS® holder

for motor cable GA804, GA805, GA806
(set consisting of 2 pieces)

56 mm



GB076R

ECCOS® holder

for perforator driver GA822

56 mm



GB077R

ECCOS® holder

for intra motor GA824

56 mm



GB082R

ECCOS® holder

for 1 handpiece

56 mm




GB083R

ECCOS® holder

for 3 handpieces

56 mm



GB085R

ECCOS® holder

for craniotome GA849

56 mm



GB086R

ECCOS® holder

for micro sagittal saw GA836

56 mm



GB087R

ECCOS® holder

for reciprocating saws GA832, GA837

56 mm



GB088R

ECCOS® holder










for saws GA831, GA833

56 mm

















AESCULAP® Mechanical Processing

ECCOS®
ELAN 4 electro

























GB482R	GB715R	GB716R
ECCOS® holder for ELAN 4 electro drill GA844  106 mm	ECCOS® holder for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R  106 mm	ECCOS® holder for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R  106 mm
 	 	 

ECCOS® & rinsing devices
ELAN 4 electro




GB084R	GB723R
ECCOS® holder for MIS handpiece GA860  56 mm	ECCOS® holder for MIS handpiece shafts GB920R – GB927R  56 mm
 	 
GB692R	GB679R
Rinsing device for app lied parts, incl. ceramic filter GB693  (not suitable for autoclaving)	Rinsing device for 4 MIS shafts, incl. ceramic filter GB693  (not suitable for autoclaving)
 	 
GB698R	GB699R
Rinsing adapter for applied parts (for manual cleaning)	Rinsing adapter for MIS shafts (for manual cleaning)
	




AESCULAP® Mechanical Processing

ECCOS®
ELAN 4 accessories

<div>GB718R</div> <div>ECCOS® holder for 12 1- and 2-ring burrs</div> <div> 56 mm</div> <div></div> <div></div>	<div>GB719R</div> <div>ECCOS® holder for 3 craniotome attachments</div> <div> 56 mm</div> <div></div> <div></div>	<div>GB720R</div> <div>ECCOS® universal holder for saw blades</div> <div> 56 mm</div> <div></div> <div></div>
<div>GB729R</div> <div>ECCOS® holder for 2 reciprocating saw shafts</div> <div> 56 mm</div> <div></div> <div></div>	<div>GB730R</div> <div>ECCOS® holder for 2 transversal saw shafts</div> <div> 56 mm</div> <div></div> <div></div>	<div>GB072R</div> <div>ECCOS® cable/hose pos. guide for cable and hoses (set of 4 pieces)</div> <div> 56 mm</div> <div></div> <div></div>
<div>GB727R</div> <div>ECCOS® holder for handpiece holding-sleeves GB861P, GB862P, GB863P</div> <div> 56 mm</div> <div></div> <div></div>	<div>GB572R</div> <div>ECCOS® fixation for GB300R - GB304R</div> <div> 56 mm</div> <div></div> <div></div>	

ECCOS® sets
Acculan 4



















<div>GB257R</div> <div>ECCOS® set for Acculan 4 small drill GA344</div> <div> 106 mm</div> <div>consisting parts: GB482R GB487R GB488R GB715R GB716R 2x GB690R JG342 JG344</div>
<div>GB243R</div> <div>ECCOS® set for 2 Acculan 4 machines</div> <div> 106 mm</div> <div>consisting parts: GB481R GB487R GB488R GB495R GB496R GB497R 2x GB690R JG342 JG344</div>
<div>GB256R</div> <div>ECCOS® set for Acculan 4 dermatome</div> <div> 106 mm</div> <div>consisting parts: GB487R GB488R GB498R JJ464R</div>


















AESCULAP® Mechanical Processing

ECCOS®














Acculan 4

GB482R	GB715R	GB716R
ECCOS® holder for Acculan 4 small drill GA344  106 mm	ECCOS® holder for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R  106 mm	ECCOS® holder for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R  106 mm
		
		
GB495R	GB496R	GB497R
ECCOS® holder for GA330, GA331, GA334  106 mm	ECCOS® holder for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R  106 mm	ECCOS® holder for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R  106 mm
		
		







GB487R	GB488R	GB489R
ECCOS® holder for GA675  56 mm	ECCOS® holder for GA678  76 mm	ECCOS® holder for GA679  106 mm
		
		
GB498R	GB686R	
ECCOS® holder for GA336, GA337, GA340, GA341  106 mm	ECCOS® holder for GB436R  56 mm	
		
		

AESCULAP® Mechanical Processing

ECCOS® holder
Acculan 4











<div><div>GB485R</div><div>ECCOS® holder for GA875</div><div> 56 mm</div></div> <div></div> <div></div>	<div><div>GB486R</div><div>ECCOS® holder for GA878</div><div> 56 mm</div></div> <div></div> <div></div>	<div><div>GB688R</div><div>ECCOS® holder for BA720R</div><div> 106 mm</div></div> <div></div> <div></div>
<div><div>GB481R</div><div>ECCOS® holder for saw blades and other accessories</div><div> 56 mm</div></div> <div></div> <div></div>	<div><div>JF148R</div><div>Mini tray with lid</div><div>outer size: 145 x 44 x 25 mm inner size: 140 x 39 x 18 mm for storage of Acculan 4 Dermatome rod and flaps during automatic cleaning process</div></div> <div></div>	

ECCOS® holder
Universal tools















<div>GB459R</div> <div>ECCOS® & Aesculap InOS® mounting kit</div> <div>consisting of: socket wrench AF7 and hex-wrench</div> <div></div>	<div>GB690R</div> <div>ECCOS® holder</div> <div>5 mm for drilling wire guide GA672207, GA344244</div> <div> </div>	<div>JG303R</div> <div>Replacement screws-set</div> <div>screw nut size AF7</div> <div></div>
<div>JG342</div> <div>Silicone clamp</div> <div>8.5 mm screw JG303R to be ordered separately</div> <div></div>	<div>JG344</div> <div>Silicone clamp</div> <div>12 mm screw JG303R to be ordered separately</div> <div></div>	

AESCULAP® Power Systems

Tool attachments overview

 1:1	AO large	 1:1	AO small
 2:1	Dental	 1:1	DIN
 1:1	Harris	 1:1	Hudson
 1:1	Hexagon	 1:1	Synthes DHS/DCS
 1:1	Trinkle	 1:1	Zimmer

Signs & symbols

	Latex-free All products listed in this catalogue are latex-free
	Article code ECCOS® fixation
	Forward/reverse drive
	Forward drive
	Oscillation frequency
	Oscillation angle
	Stroke
	Cannulation
	Cable length
	Weight
	Packaging unit
	Single use
	Required height of tray
	Sterile


AESCULAP® Power Systems

Brochures

Index

Main Catalogues


084002
Power Systems



017599
Burrs & Blades



CLICK HERE
for Electronic Instructions for Use (eIFU)




ELAN 4

071602
ELAN 4

071511
Mechanical Cleaning






012002 Manual Cleaning





Acculan 4


001802
Acculan 4


001902
Acculan 4 Dermatome002002
Acculan 4 Electro


002102
Acculan 4 Sternum Saws



Art. No.	
BA720R	32
GA031R	42
GA062R	42
GA258SU	23
GA259SU	23
GA261SU	23
GA262SU	23
GA263SU	23
GA264SU	23
GA265SU	23
GA266SU	23
GA267SU	23
GA268SU	23
GA269SU	23
GA320	38
GA330	26
GA331	30
GA334	31
GA336	30
GA337	30
GA340	32
GA341	32
GA344	34
GA346	38
GA395SU	23
GA675	38
GA676	38
GA678	38, 41
GA679	38
GA800	7
GA804	8
GA805	8
GA806	7
GA808	7
GA810	7
GA822	8
GA824	8
GA831	15
GA832	15
GA833	15
GA836	15
GA837	15
GA844	17
GA849	9
GA860	13
GA861	10
GA862	10
GA863	10
GA864	10
GA865	11

Art. No.	
GA866	11
GA867	11
GA868	11
GA869	11
GA875	41
GA876	41
GA877	41
GA878	41
GB072R	50
GB073R	47
GB076R	47
GB077R	47
GB082R	47
GB083R	47
GB084R	49
GB085R	47
GB086R	47
GB087R	47
GB088R	47
GB176R	42
GB184R	42
GB243R	51
GB256R	51
GB257R	51
GB261R	46
GB262R	46
GB414R	42
GB436R	31
GB459R	55
GB481R	54
GB482R	48, 52
GB485R	54
GB486R	54
GB487R	53
GB488R	53
GB489R	53
GB495R	52
GB496R	52
GB497R	52
GB498R	53
GB572R	50
GB600	44
GB600840	44
GB600860	44
GB600870	44
GB600880	44
GB620R	26
GB621R	26
GB623R	26
GB628R	42

Art. No.	
GB639R	27
GB641R	29
GB645R	27
GB657R	28
GB660R	29
GB663R	27
GB664R	27
GB665R	27
GB667R	28
GB668R	28
GB669R	28
GB670R	28
GB679R	49
GB686R	53
GB688R	54
GB690R	55
GB692R	49
GB698R	49
GB699R	49
GB715R	48, 52
GB716R	48, 52
GB718R	50
GB719R	50
GB720R	50
GB723R	49
GB727R	50
GB729R	50
GB730R	50
GB796SU	23
GB797SU	23
GB861P	22
GB862P	22
GB863P	22
GB866R	22
GB870R	17, 34
GB871R	17, 34
GB872R	17, 34
GB878R	19, 36
GB881R	18, 35
GB882R	18, 35
GB884R	18, 35
GB886R	18, 35
GB887R	18, 35
GB891R	19, 36
GB892R	19, 36
GB894R	19, 36
GB896R	20, 37
GB897R	20, 37
GB920R	13
GB921R	13

Art. No.	
GB922R	13
GB925R	13
GB926R	13
GB927R	13
GB941R	9
GB942R	9
GB943R	9
GB945R	9
GB947R	9
JF148R	54
JF212R	46
JG303R	55
JG342	55
JG344	46, 55
TE730	22, 42
TE734	22, 42
TE735	22, 42
TE780	22, 42
XF457R	42

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

Distributed by

B. Braun Medical B.V. | +31 (0)412 67 24 11 | customercare.nl@bbraun.com | www.bbraun.nl

The main product trademark "AESCULAP" and the product trademarks "Aesculap InOS", "ECCOS" and "Sterilit" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.