

Easypump® II

Over- en ondervullen

Productspecificaties

Easypump® II ST										
Short Term – variant	100 – 0.5	250 – 0.5	50 – 1	100 – 1	250 – 1	250 – 1.5	400 – 2	500 – 2	100 – 2	400 – 4
Artikelnummer	4540040-07	4540042-07	4540044-07	4540046-07	4540048-07	4540050-07	4540052-07	4540054-07	4540056-07	4540058-07
Nominale toedieningssnelheid (ml/u)	200	500	50	100	250	175	200	250	50	100
Nominaal vulvolume (ml)	100	250	50	100	250	250	400	500	100	400
Minimaal vulvolume (ml)	50	125	30	50	135	135	240	240	50	240
Maximaal vulvolume (ml)	125	295	65	125	295	295	560	560	125	560
Maximaal resterend volume (ml)	≤3	≤8	≤2	≤3	≤8	≤8	≤10	≤10	≤3	≤10

Bewaren

Tijd om kamertemperatuur te bereiken (circa +23°C ±2°C)										
Aantal uren uit de koelkast (circa +2°C tot +8°C)	4	6	4	4	6	6	8	8	4	8
Aantal uren uit de vriezer (circa -18°C)	12	18	12	12	18	18	18	18	12	18

Data over- en ondervullen

Geschatte toedieningstijd										
Uren	Vulvolume (ml)									
00:15 uur	50	125	-	-	-	-	-	-	-	-
00:30 uur	100	250	-	50	-	-	-	-	-	-
00:45 uur	-	-	38	75	188	-	-	-	-	-
01:00 uur	-	-	50	100	250	175	-	250	50	-
01:15 uur	-	-	63	125	-	219	-	313	63	-
01:30 uur	-	-	-	-	-	263	300	375	75	-
01:45 uur	-	-	-	-	-	-	350	438	88	-
02:00 uur	-	-	-	-	-	-	400	500	100	-
02:15 uur	-	-	-	-	-	-	450	-	113	-
02:30 uur	-	-	-	-	-	-	500	-	125	250
02:45 uur	-	-	-	-	-	-	550	-	-	275
03:00 uur	-	-	-	-	-	-	-	-	-	300
03:30 uur	-	-	-	-	-	-	-	-	-	350
04:00 uur	-	-	-	-	-	-	-	-	-	400
04:30 uur	-	-	-	-	-	-	-	-	-	450
05:00 uur	-	-	-	-	-	-	-	-	-	500
05:30 uur	-	-	-	-	-	-	-	-	-	550



Let op: om het nominale vulvolume toe te dienen, moet je het maximaal resterend volume toevoegen.
Voorbeeld ST 100-0.5:
100 ml + 3 ml = 103 ml (nominaal vulvolume + maximaal resterend volume = totale vulvolume).

Productspecificaties

Easypump® II LT																				
Long Term – variant	60 – 12	500 – 12.5	80 – 16	125 – 25	270 – 27	60 – 30	120 – 30	400 – 40	100 – 50	270 – 54	120 – 60	400 – 80	100 – 67	270 – 68	400 – 100	65 – 130	270 – 135	300 – 150	100 – 200	270 – 270
Artikelnummer	4540002-07	4540003-07	4540004-07	4540006-07	4540008-07	4540010-07	4540012-07	4540014-07	4540016-07	4540018-07	4540020-07	4540022-07	4540024-07	4540026-07	4540028-07	4540030-07	4540032-07	4540034-07	4540036-07	4540038-07
Nominale toedieningssnelheid (ml/u)	5	40	5	5	10	2	4	10	2	5	2	5	1.5	4	4	0.5	2	2	0.5	1
Nominaal vulvolume (ml)	60	500	80	125	270	60	120	400	100	270	120	400	100	270	400	65	270	300	100	270
Minimaal vulvolume (ml)	30	240	50	60	120	30	60	240	50	135	60	240	50	135	240	30	135	150	50	135
Maximaal vulvolume (ml)	65	560	125	125	295	65	125	560	125	295	125	560	125	295	560	65	295	340	125	295
Maximaal resterend volume (ml)	≤2	≤10	≤3	≤3	≤8	≤2	≤3	≤10	≤3	≤8	≤3	≤10	≤3	≤8	≤10	≤2	≤8	≤9	≤3	≤8

Bewaren

Tijd om kamertemperatuur te bereiken (circa +23°C ±2°C)																				
Aantal uren uit de koelkast (circa +2°C tot +8°C)	4	8	4	5	6	4	5	8	4	6	5	8	4	6	8	4	6	6	4	6
Aantal uren uit de vriezer (circa -18°C)	12	18	12	18	18	12	18	18	12	18	18	18	12	18	18	12	18	18	12	18

Data over- en ondervullen

Geschatte toedieningstijd																				
Uren	Dag(en)	Vulvolume (ml)																		
6 uur	-	-	240	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8 uur	-	40	320	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10 uur	-	50	400	50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12 uur	-	60	480	60	60	120	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18 uur	-	-	-	90	90	180	36	72	-	-	-	-	-	-	-	-	-	-	-	-
24 uur	1 dag	-	-	120	120	240	48	96	240	-	-	-	-	-	-	-	-	-	-	-
30 uur	-	-	-	-	-	-	60	120	300	60	150	60	-	-	-	-	-	-	-	-
48 uur	2 dagen	-	-	-	-	-	-	-	480	96	240	96	240	72	-	-	-	-	-	-
60 uur	2.5 dagen	-	-	-	-	-	-	-	-	120	-	120	300	90	240	240	30	-	-	-
72 uur	3 dagen	-	-	-	-	-	-	-	-	-	-	-	360	108	288	288	36	144	-	-
96 uur	4 dagen	-	-	-	-	-	-	-	-	-	-	-	480	-	-	384	48	192	192	-
120 uur	5 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	480	60	240	240	60
144 uur	6 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	288	288	72
168 uur	7 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	336	84
192 uur	8 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	96
216 uur	9 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	108
240 uur	10 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	120
264 uur	11 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
288 uur	12 dagen	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

Let op: om het nominale vulvolume toe te dienen, moet je het maximaal resterend volume toevoegen.
Voorbeeld LT 60-12: 60 ml + 2 ml = 62 ml (nominaal vulvolume + maximaal resterend volume = totale vulvolume).