



## **Venostrip**®

### Disposable vein extirpation set

#### Introduction

Venous diseases in legs belong to the most common found diseases. They exceed by far all other vascular diseases of the leg. In developed countries approximately 10 – 20 % of the population suffer from varicose veins and appr. 0.5 % from ulcers of the leg.<sup>1, 2, 3</sup>
For successful treatment of this frequent disease the easy handling of the instruments is most important.
For this purpose the Venostrip® is a vein extirpation set offers the following advantages.

#### Product description

Venostrip®, a disposable vein stripper according to Gruß, consists of a plastic-coated guidewire, a handle and three stripping heads in 9, 12 and 15 mm diameter and one screw cap to fasten the stripping heads.

- Yang L, Wang XP, Su WJ, Zhang Y, Wang Y. Randomized clinical trial of endovenous microwave ablation combined with high ligation versus conventional surgery for varicose veins. Eur J Vasc Endovasc Surg 2013; 46(4):473–9.
- Velička L, Kievišas M, Keturakis V. A comparison of endovenous laser ablation and conventional surgery in patients with varicose veins of the lower limbs. Acta Angiologica 2016; 21(4):107–15.
- <sup>3</sup> Ran F, Shi Y, Qiao T, Shang T, Liu Z, Liu C-J. Comparison of Foam Sclerotherapy Alone or Combined With Stripping of the Great Saphenous Vein for Treating Varicose Veins. Dermatol Surg 2017; 43(4):541–7.





### Advantages

- Polyfil, flexible, plastic-coated guidewire:
  - the wire offers flexible handling properties
  - no curling of the wire leads to precise guidance
  - for external and internal stripping
- Sharp-edged stripping heads:
  - diameter 9, 12, and 15 mm
  - improved cosmetics results compared to the preoperative condition
- Single use product:
  - medical safety aspect for patients and users

### Intended use & Indication

Extirpation of the great and small saphenous vein or parva vein in primary varicosis, chronic venous insufficiency and in the case of serious cosmetic varicosity. Venostrip® can be used for both exostripping and endostripping.

### Notes



# Venostrip<sup>®</sup>

### Disposable vein extirpation set



### Product range

Article No.	Contents
5523184	Venostrip® Complete set, sterile Carton of 10 pieces
5523176	Venostrip® Stripper wire, sterile Carton of 10 pieces





#### Manufacturer

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany CE  $\,$  0123

### Distributor

B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany Phone +49 30 568207-300 | Fax +49 30 568207-130 | www.bbraun.com

The product trademark "Venostrip" is a registered trademark of Aesculap AG. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.