

# ECOMIX – CONCENTRATE PORTFOLIO

## ECOCart REVOLUTION Portfolio Acid Concentrates, with Acetate 1+44

Article	Product name	K <sup>+</sup> mmol/l	Na <sup>+</sup> mmol/l	Ca <sup>++</sup> mmol/l	Mg <sup>++</sup> mmol/l	Cl <sup>-</sup> mmol/l	HCO <sub>3</sub> <sup>-</sup> mmol/l	Acetate mmol/l	Glc. g/l	Theor.Osmol. mOsm/l
IM3655	DSK-3252-GW-622-45-Rev	1	138	1.25	0.5	107.5	32	3	1	288.8
IM3726	DSK-3262-GW-622-45-Rev	1	138	1.5	0.5	108	32	3	1	289.55
IM3730	DSK-3272-GW-622-45-Rev	1	138	1.75	0.5	108.5	32	3	1	290.3
IM3452	DSK-3452-GW-622-45-Rev	2	138	1.25	0.5	108.5	32	3	1	290.8
IM3453	DSK-3462-GW-622-45-Rev	2	138	1.5	0.5	109	32	3	1	291.55
IM3952	DSK-3472-GW-622-45-Rev	2	138	1.75	0.5	109.5	32	3	1	292.3
IM3635	DSK-3652-GW-622-45-Rev	3	138	1.25	0.5	109.5	32	3	1	292.8
IM3750	DSK-3662-GW-622-45-Rev	3	138	1.5	0.5	110	32	3	1	293.55
IM3957	DSK-3672-GW-622-45-Rev	3	138	1.75	0.5	110.5	32	3	1	294.3
IM3962	DSK-3852-GW-622-45-Rev	4	138	1.25	0.5	110.5	32	3	1	294.8
IM3937	DSK-3862-GW-622-45-Rev	4	138	1.5	0.5	111	32	3	1	295.55
IM3968	DSK-3872-GW-622-45-Rev	4	138	1.75	0.5	111.5	32	3	1	296.3
IM4005	DSK-3253-GW-622-45-Rev	1	138	1.25	0.75	108	32	3	1	289.55
IM4009	DSK-3263-GW-622-45-Rev	1	138	1.5	0.75	108.5	32	3	1	290.3
IM4013	DSK-3273-GW-622-45-Rev	1	138	1.75	0.75	109	32	3	1	291.05
IM4017	DSK-3453-GW-622-45-Rev	2	138	1.25	0.75	109	32	3	1	291.55
IM4020	DSK-3463-GW-622-45-Rev	2	138	1.5	0.75	109.5	32	3	1	292.3
IM4024	DSK-3473-GW-622-45-Rev	2	138	1.75	0.75	110	32	3	1	293.05
IM3748	DSK-3653-GW-622-45-Rev	3	138	1.25	0.75	110	32	3	1	293.55
IM4028	DSK-3663-GW-622-45-Rev	3	138	1.5	0.75	110.5	32	3	1	294.3
IM4032	DSK-3673-GW-622-45-Rev	3	138	1.75	0.75	111	32	3	1	295.05
IM3995	DSK-3853-GW-622-45-Rev	4	138	1.25	0.75	111	32	3	1	295.55
IM4038	DSK-3863-GW-622-45-Rev	4	138	1.5	0.75	111.5	32	3	1	296.3
IM4042	DSK-3873-GW-622-45-Rev	4	138	1.75	0.75	112	32	3	1	297.05

## ECOCart REVOLUTION Portfolio Acid Concentrates, with Citrate, Acetate-free 1+44

Article	Product name	K <sup>+</sup> mmol/l	Na <sup>+</sup> mmol/l	Ca <sup>++</sup> mmol/l	Mg <sup>++</sup> mmol/l	Cl <sup>-</sup> mmol/l	HCO <sub>3</sub> <sup>-</sup> mmol/l	Citrate mmol/l	Glc. g/l	Theor.Osmol. mOsm/l
IM3768	DSK-3462-GW-Cit-622-45-Rev	2	138	1.5	0.5	109	32	1	1	289.55
IM3820	DSK-3472-GW-Cit-622-45-Rev	2	138	1.75	0.5	109,5	32	1	1	290.3
IM3827	DSK-3662-GW-Cit-622-45-Rev	3	138	1.5	0.5	110	32	1	1	291.55
IM3822	DSK-3672-GW-Cit-622-45-Rev	3	138	1.75	0.5	110,5	32	1	1	292.3
IM3826	DSK-3862-GW-Cit-622-45-Rev	4	138	1.5	0.5	111	32	1	1	293.55